

THE BOBO WOBBLY SUPPORT CUSHION

User guide


BoBo Wobbly **Support Cushion** allows you to adjust the difficulty level of your training while you use the balance board. The cushion makes your training easier by preventing the device from spinning and limiting the tilting angle.



- 1 Inflate the support cushion. The more you inflate it, the more stable and easy to use the balance board will be.



- 2 Place the support cushion on a flat, level floor and clear the training area of any hazards.



- 3 Place the BoBo Wobbly Balance Board on top and exactly in the center of the support cushion.

To prevent injury, please read **the full BoBo Wobbly Balance Board user guide** carefully before you use it



*Maximum weight - 150 kg / 330 lb

 **bobo**
Balance Better.