



T37 • 4655004198

Montage complet sur Lit Linet Eleganza Smart
Complete set on Linet Eleganza Smart bed

Adaptation spéciale pour lit :
Special Adaptation for bed:



Dimensions en mm.
Sizes in mm.



T374655004198Montage completComplete set

Le cadre est composé de :

The frame consists of these parts:

PRT T374655004213	Portique Tête de lit Linet Eleganza Smart (x1)	Head porch bar for Linet Eleganza Smart (x1)
PRP T374655004221	Portique Pied de lit Eleganza Smart (x1)	Foot porch bar for Eleganza Smart (x1)
T1804655002978	Barre octogonale de 180cm (x1)	Plain bar of 180cm (71") (x1)
DME Equi4655003059	Double mâchoire (x1)	Double clamp (x1)
T 40 Equi4655003067	Tube à mâchoire fixe de 40cm (x1)	Fixed clamp bar of 40cm (16") (x1)
P 7534655000948	Etrier de traction (x1)	Handgrip (x1)
PPO4655000617	Porte poulie ouverte (x7)	Open pulley with clamp (x7)
SP 014635001304	Sphère de protection (x2)	Wall protection cushion (x2)

Construit à partir d'un matériau qui allie la résistance et la légèreté, le Système de Traction KINETEC® est d'un emploi facile pour tout le personnel soignant.

Designed from a material which combines resistance and lightness, the KINETEC® TRACTION System is easy to use for all medical staff.

Kinetec SAS garantit son Système de Traction, 2 ans contre tous vices de fabrication, à partir de la date d'achat par le consommateur.

Kinetec SAS guarantees its Traction System 2 years against all defects of manufacture, starting from the date of purchase by the consumer.

Dispositif médical conforme aux exigences essentielles de la directive 93/42/CEE relative aux dispositifs médicaux de classe I.

Medical device conforms to the essential requirements of the directive 93/42/CEE relating to the medical devices of class I.

Demandez notre catalogue complet d'accessoires et de kits de traction (poids, corde, attelle de positionnement, etc...)

Ask for our complete catalogue of accessories and traction kits (weights, cord, positioning splints, and more...)